Coral Clinical Systems

Plot No. M – 46, Phase III B, Verna Ind. Estate,
Verna, Goa – 403 722, INDIA
Tel.: (0832) 6680121, Fax: (0832) 2887028

Corporate Office : "Gitanjali", Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz,
Bambolim Complex Post Office, Goa – 403 202. INDIA.
Tel.: (0832) 2458546 – 51, Fax: (0832) 2458544

Material Safety Data Sheet

1. Product Identification

Product Name : **Calkine CK MB (NAC Act.) Kit (Immunoinhibition/Mod. IFCC Method)**

Catalog Number : CKB 110 / CKB 111

2. Composition / Information on Hazardous Ingredients

Chemical Name	CAS #	% W/V	Exposure Limits in Air				
			ACGIH		OSHA		OTHER
			TLV	STEL	PEL	STEL	
Reagent 1: Enzyme Reagent							
Sodium Azide	26628-22-8	0.1	N / A	0.29 mg/m^3	N / A	N / A	NIOSH 0.3 mg/m^3 C (skin)
Reagent 2: Starter Reagent							
Sodium Azide	26628-22-8	0.1	N / A	0.29 mg/m^3	N / A	N / A	NIOSH 0.3 mg/m^3 C (skin)
Reagent 3: N/A							
Reagent 4: N/A							
Reagent 5: N/A							

3. Hazard Identification

Primary Routes of Entry:
Inhalation, ingestion, skin and / or eye contact.

Inhalation:
Inhalation of Vapours, mists, or sprays of these components may irritate the nose, throat, and lungs. Symptoms are generally alleviated upon breathing fresh air.

Ingestion:
Though not a likely route of occupational exposure, ingestion of this product, especially in large quantities, May cause gastric distress. Symptoms may include nausea, vomiting, or diarrhea.

Skin Contact:
Sodium Azide can be absorbed through skin. Skin contact may cause reddening, discomfort, or mild irritation. Irritation should be alleviated upon rising.

Plot No. M – 46, Phase III B, Verna Ind. Estate,
Verna, Goa – 403 722, INDIA
Tel.: (0832) 6680121, Fax: (0832) 2887028

Eye Contact:
May cause irritation.
Chronic Exposure:
N/A

Medical Conditions Aggravated by Exposure:
N/A

Health Effects:
No adverse health effects should occur from routine use of this material in the manner specified by area work instructions.

4. First Aid Measures

Inhalation:
If breathing becomes difficult, remove victim to fresh air. Seek medical attention immediately.

Ingestion:
Do not induce vomiting. Get medical attention if there has been ingestion of this product.

Skin Contact:
Avoid skin contact. If skin contact occurs, remove contaminated clothing and wash exposed skin with water for atleast 15 minutes. Get medical attention immediately.

Eye Contact:
Immediately flush eye(s) with large volume of water for atleast 15 minutes, occasionally lifting the lower lids. Get medical attention immediately.

5. Fire Fighting Measures

Flash Point (Method used): N/A **Flammable Limits – LEL:** N/A **UEL:** N/A

Extinguishing Media:
Use fire extinguishing media appropriate for site conditions.

Special Fire Procedures:
Structural firefighting gear and self-contained breathing apparatus will provide adequate protection if this product is in a fire area.

Unusual Fire and Explosion Hazards:
Sodium azide can react with copper, lead, brass, or solder in plumbing to form explosive compound of lead azide and copper azide. Sodium azide can react with acids to form explosive hydrogen azide.

6. Accidental Release Measures

Steps to be taken in case material is Released or Spilled:
PPE should be level D: lab gloves, chemical resistant apron, boots and splash goggles. Use an absorbent material to contain / pick up the spilled solution. Place all spill residue into a suitable container, seal, label and hold for disposal.

7. Handling and Storage

Refer to packet insert for additional information on handling and storage procedures.

8. Exposure Controls and Personal Protection

Ventilation Data:
A system of local and / or general exhaust is recommended to keep employee exposures as low as possible.

Respiratory Protection:
Respiratory protection is not required under normal use of this product. If respiratory protection is needed, follow OSHA respirator regulations (29CFR1910.134) and, if necessary, wear a NIOSH approved respirator. Select respirator based on its suitability to provide worker protection for given working conditions, level of airborne concentration, and presence of sufficient oxygen.

Plot No. M – 46, Phase III B, Verna Ind. Estate,
Verna, Goa – 403 722, INDIA
Tel.: (0832) 6680121, Fax: (0832) 2887028

Protective Gloves:
Wear appropriate protective gloves to prevent skin contact. Replace torn or punctured gloves promptly.

Other Protective Equipment:
Wear appropriate eye protection to prevent eye contact. Wear appropriate body protection to prevent skin contact.

Other Engineering Controls:
Eye wash stations and deluge showers.

Work Practices:
Good laboratory technique should be used when handling this product. Observe appropriate chemical hygiene. Avoid contact with eyes or skin. Do not place in mouth.

Hygienic Practices:
Do not eat, drink, or smoke while working with product. Upon completion of work activities involving this product, wash hands thoroughly with soap and water.

9. Physical And Chemical Properties

For All Components Unless Otherwise Indicated

Relative Vapour density(air = 1)	:	N/A	Evaporation rate(nBuAc = 1):		N/A
Specific Gravity (water = 1)	:	N/A	Freezing / Melting Point	:	N/A
Solubility in Water	:	Soluble	Boiling Point	:	N/A
Vapour Pressure, mm Hg @ 20ºC:		N/A	pH	:	N/A

Odour and Appearance Information
Reagent 1: pH=6.05, Clear, Colourless liquid

Reagent 2: pH=9.2, Clear, Colourless liquid

Reagent 3: N/A

Reagent 4: N/A.

Reagent 5: N/A.

10. Stability and Reactivity

Incompatibility (Materials to Avoid):
Strong bases, strong acids, and water reactive materials.

Hazardous Decomposition Products:
Thermal decomposition may emit toxic products.

Will Hazardous Polymerization Occur?
Hazardous polymerization will not occur.

Conditions to Avoid / Polymerization: N/A
Is the Product Stable?
Yes, under normal handling and storage conditions.
Conditions to Avoid/stability
Stable Solution. Avoid acidification of solution, which may generate hydrogen cyanide gas.

11. Toxicological Information

Toxicity Data:
Sodium Azide (Undiluted): LD50 (rat and mouse, Oral)=27 mg/kg; LD50 (rabbit, skin)=20 mg/kg.
Reproductive effects:
N/A.

Plot No. M – 46, Phase III B, Verna Ind. Estate,
Verna, Goa – 403 722, INDIA
Tel.: (0832) 6680121, Fax: (0832) 2887028

Target organ Effects:
Eyes, Skin, central nervous systems, cardiovascular systems, and digestive System.

Carcinogenicity: No

CHEMICAL NAME	CAS #	% W/V	NTP Carcinogen		IARC	OSHA
			Known	Anticipated		
N/A.						

12. Ecological Information

Environmental Fate / Stability:
N/A

Effect of Material on plants or animals:
N/A

Effect of Chemical on Aquatic Life:
N/A

13. Disposal Considerations

EPA Waste Number and Proper Waste Disposal Method:
Please consult local, state and federal regulations for additional guidance on disposal.

14. Transportation Information

Is this Material Hazardous? Not regulated under transportation regulations.

Proper Shipping Name : N/A		
Hazard Class Number : N/A	Packing Group: N/A	UN Number: N/A

15. Regulatory Information

NA.

16. Other Information

NA => NOT APPLICABLE or NO INFORMATION