
ALPHA AMYLASE

Optically Clean Sterile Filled Assay System



ALPHA AMYLASE



CORAL ALPHA AMYLASE KIT

Amylase reactions during assay are very sensitive to pH.

Liquid stable preparations ensure that pH is buffered to the optimum

Substrate is GAL G2 CNP

Sensitive to both Salivary and Pancreatic Amylase thus avoiding underestimation of Total Amylase

Reagent contains optimum Ca2+ ions

Maintains functional integrity of amylase in sample

Unaffected by high Triglycerides or Bilirubin

Can be used even with lipemic and/or icteric samples

Manufactured by :

Coral Clinical Systems

TULIP GROUP OF COMPANIES
Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403202, INDIA.
Tel. : + 91 832 2458546-50 **Fax :** + 91 832 2458544 **E-mail :** sales@tulipgroup.com **Website :** www.tulipgroup.com