
CK NAC & CK MB



- Thiol Activator to reactivate CK in specimen
- Optimised EDTA
- Optimised blend of Adenylate Kinase inhibitors
- Purified Fab fragment-based Anti CK-M antibody
- Optically clean liquid reagent assay system

CK NAC & CK MB





Features	Benefits
■ Optimum Conc. of THIOL ACTIVATOR	For the maintenance and reactivation of CK activity in stored samples
■ OPTIMIZED EDTA	To avoid interference of cations
■ OPTIMIZED BLEND of Adenylate kinase inhibitors	To avoid increased or false values of CK
■ PURIFIED FAB FRAGMENT-BASED Anti CK – M antibody	To avoid Fc fragment interference to deliver accurate results
■ OPTICALLY CLEAN liquid reagent assay system	No interference from suspended micro particles

CK Nac & CK MB....
The ideal Cardiac Markers

Manufactured by :

Coral Clinical Systems

orchid · Microxpress · Coral Clinical Systems · BioShields · Viola · molbio

TULIP GROUP OF COMPANIES

Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403202, INDIA.

Tel. : + 91 832 2458546-50 **Fax :** + 91 832 2458544 **E-mail :** sales@tulipgroup.com **Website :** www.tulipgroup.com