
Discover the way of Healthy Living




Check Your Triglycerides Levels at Regular Interval

TRIGLYCERIDE





TRIGLYCERIDES HIGH LINEARITY



TRIGLYCERIDES - LINEARITY

Conc.	0 min	10 min
200 mg/dl	0	0.240
400 mg/dl	0	0.478
600 mg/dl	0	0.717
800 mg/dl	0	0.948
1000 mg/dl	0	1.175

CORAL TRIGLYCERIDES KIT

LIQUID STABLE ready to use kit

ENZYME and chromogen provided separately, Triglyceride SR (single reagent) kit also available

RAPID and stable end points

LARGER SHELF LIFE and lower blanks

HIGH linearity (1000 mg/dl)

CONVENIENT pack sizes

Manufactured by.

Coral Clinical Systems

TULIP GROUP OF COMPANIES

Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403202, INDIA.

Tel. : + 91 832 2458546-50 **Fax :** + 91 832 2458544 **E-mail :** sales@tulipgroup.com **Website :** www.tulipgroup.com