# CORALAB® Touch
CLINICAL CHEMISTRY ANALYZER

| TECHNICAL SPECIFICATIONS | |
|---|---|
| **Operating Modes** | Absorbance, End-point, Fixed Time, Kinetic, Multipoints, Sample blank |
| **Items** | 180 Test Parameters |
| **Result Storage** | Upto 6000 Test Results |
| **Reaction Volume** | 400 µl |
| **Flow Cell** | 32 µl |
| **Temperature Control** | 37°C ± 0.2°C |
| **Optical System** | 340, 405, 492, 510, 546, 578, 620 nms |
| **Light Source** | 6V, 10W Halogen Lamp |
| **Photometric Range** | 0 – 3.0 A |
| **Software** | Touch Screen, English Version |
| **Display** | 7 inch Touch Screen |
| **Printer** | Built – in Thermal Printer |
| **Serial Output** | RS – 232 Standard |
| **Power Supply** | 220V/110V ± 10% 50 Hz – 60 Hz |
| **Dimension** | 350mm (W) x 390mm (L) x 184mm (H) |
| **Weight** | 6 kgs |

** Technical features are subject to change without prior notice.*

Product of:

Coral | Clinical Systems

A Division of

**TULIP DIAGNOSTICS (P) LTD.**

CIN: U85195GA1988PTC005717

Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403202, INDIA.
**Tel.:** +91 832 2458546-50 **Fax :** +91 832 2458544 **E-mail :** sales@tulipgroup.com **Website :** www.tulipgroup.com

CLTLT10-19/VER-01

Coral | Clinical Systems

CORALAB® Touch
CLINICAL CHEMISTRY ANALYZER

"Empowered by Innovation"






**Be in Touch with better performance...**


Size : 420mm W x 297mm H

# CORALAB® Touch


## CLINICAL CHEMISTRY ANALYZER

*Coralab Touch is a semi-automatic analyzer with a 7 inch touchscreen display, 20 thermostat positions programmable from 25°C to 37°C. It stores upto 180 test parameters and 6000 test results. It is a fully open system. Optics and automatic filter selection assures repeatable measurements. Results are displayed and printed out directly in the required measuring unit.*

### Programming Mode





### Main Menu





Upto 180 test parameters




### Levey Jennings Graph





### Real Time Graph




Real time plot for kinetic & Fixed Time


### Temperature Selection




Adaptive to various Chemistries




### Sample Position



20 Thermostat Positions

***Be in Touch***
***with better performance...***

Size : 420mm W x 297mm H