ChemCAL™

Chemistry Calibrator

LOT 13925051A IVD Mfd. May. 2025 Exp. Mar. 2027

BT 1500

PRODUCT	METHOD	UNITS	CONCENTRATION
ALBUMIN	BROMOCRESOL GREEN	g/dl	4.67
ALKALINE PHOSPHATASE	IFCC (AMP)	U/L	452
ALT/ GPT	MOD. IFCC	U/L	119
AST/ GOT	MOD. IFCC	U/L	105
AMYLASE	CNPG3, KINETIC	U/L	228
BILIRUBIN DIRECT	MOD. JENDRASSIK-GROF	mg/dl	1.75
	DCA		3.14
BILIRUBIN TOTAL	MOD. JENDRASSIK-GROF	mg/dl	4.94
	DCA		5.39
CALCIUM	AIII	mg/dl	10.5
CHLORIDE	THIOCYANATE	mmol/l	201
CHOLESTEROL	CHOD-PAP	mg/dl	112
HDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	58
LDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	92
CREATINE KINASE (CK)	MOD. IFCC	U/L	241
CREATININE	MOD. JAFFE'S KINETIC	mg/dl	3.36
	ENZYMATIC PAP	mg/dl	3.89
ɣ- GLUTAMYL TRANSFERASE (GGT)	CARBOXY SUBSTRATE	U/L	122
GLUCOSE	GOD/POD	mg/dl	208
IRON	FERROZINE	µg/dl	131
LACTATE DEHYDROGENASE	DGKC (P-L)	U/L	451
LIPASE	COLORIMETRIC	U/L	100
MAGNESIUM	CALMAGITE	mEq/L	3.10
PHOSPHORUS	MOLYBDATE U. V	mg/dl	5.46
TOTAL PROTEIN	BIURET	g/dl	7.42
TRIGLYCERIDES	GPO/PAP	mg/dl	174
UREA	GLDH KINETIC	mg/dl	115
URIC ACID	URICASE /PAP	mg/dl	5.79

FOR IN VITRO DIAGNOSTIC USE ONLY

Manufactured By:
Coral Clinical Systems
A Division of Tulip Diagnostics (P) Ltd.
Building E, Plot No. M-46/47, Phase III B,
Verna Industrial Estate, Verna, Goa – 403 722, India.

REGD. OFFICE:
Gitanjali, Tulip Block,
Dr. Antonio Do Rego Bagh,
Alto Santacruz, Bambolim Complex P.O.,
Goa – 403 202, India.

ChemCAL™

Chemistry Calibrator

LOT 13925051A IVD Mfd. May. 2025 Exp. Mar. 2027

BT 4500

PRODUCT	METHOD	UNITS	CONCENTRATION
ALBUMIN	BROMOCRESOL GREEN	g/dl	4.17
ALKALINE PHOSPHATASE	IFCC (AMP)	U/L	444
ALT/ GPT	MOD. IFCC	U/L	120
AST/ GOT	MOD. IFCC	U/L	109
AMYLASE	CNPG3, KINETIC	U/L	221
BILIRUBIN DIRECT	MOD. JENDRASSIK-GROF	mg/dl	1.65
	DCA		3.11
BILIRUBIN TOTAL	MOD. JENDRASSIK-GROF	mg/dl	4.34
	DCA		5.96
CALCIUM	AIII	mg/dl	11.1
CHLORIDE	THIOCYANATE	mmol/l	193
CHOLESTEROL	CHOD-PAP	mg/dl	103
HDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	50
LDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	86
CREATINE KINASE (CK)	MOD. IFCC	U/L	233
CREATININE	MOD. JAFFE'S KINETIC	mg/dl	3.18
	ENZYMATIC PAP	mg/dl	3.32
ɣ- GLUTAMYL TRANSFERASE (GGT)	CARBOXY SUBSTRATE	U/L	117
GLUCOSE	GOD/POD	mg/dl	196
IRON	FERROZINE	µg/dl	127
LACTATE DEHYDROGENASE	DGKC (P-L)	U/L	446
LIPASE	COLORIMETRIC	U/L	95
MAGNESIUM	CALMAGITE	mEq/L	3.40
PHOSPHORUS	MOLYBDATE U. V	mg/dl	5.16
TOTAL PROTEIN	BIURET	g/dl	6.72
TRIGLYCERIDES	GPO/PAP	mg/dl	165
UREA	GLDH KINETIC	mg/dl	105
URIC ACID	URICASE /PAP	mg/dl	5.27

FOR IN VITRO DIAGNOSTIC USE ONLY

Manufactured By:
Coral Clinical Systems
A Division of Tulip Diagnostics (P) Ltd.
Building E, Plot No. M-46/47, Phase III B,
Verna Industrial Estate, Verna, Goa – 403 722, India.

REGD. OFFICE:
Gitanjali, Tulip Block,
Dr. Antonio Do Rego Bagh,
Alto Santacruz, Bambolim Complex P.O.,
Goa – 403 202, India.

ChemCAL™

Chemistry Calibrator

LOT 13925051A IVD Mfd. May. 2025 Exp. Mar. 2027

CORALYZER SMART / JUNIOR

PRODUCT	METHOD	UNITS	CONCENTRATION
ALBUMIN	BROMOCRESOL GREEN	g/dl	4.40
ALKALINE PHOSPHATASE	IFCC (AMP)	U/L	400
	DGKC (DEA)	U/L	438
ALT/ GPT	MOD. IFCC	U/L	115
AST/ GOT	MOD. IFCC	U/L	101
AMYLASE	CNPG3, KINETIC	U/L	206
BILIRUBIN DIRECT	MOD. JENDRASSIK-GROF	mg/dl	1.98
	DCA		2.81
BILIRUBIN TOTAL	MOD. JENDRASSIK-GROF	mg/dl	7.49
	DCA		4.97
CALCIUM	OCPC	mg/dl	10.4
	AIII	mg/dl	10.7
CHLORIDE	THIOCYANATE	mmol/l	107
CHOLESTEROL	CHOD-PAP	mg/dl	207
HDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	62
LDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	83
CREATINE KINASE (CK)	MOD. IFCC	U/L	244
CREATININE	MOD. JAFFE'S KINETIC	mg/dl	3.25
	ENZYMATIC PAP	mg/dl	4.23
ɣ- GLUTAMYL TRANSFERASE (GGT)	CARBOXY SUBSTRATE	U/L	122
GLUCOSE	GOD/POD	mg/dl	180
IRON	FERROZINE	µg/dl	140
LACTATE DEHYDROGENASE	DGKC (P-L)	U/L	494
LIPASE	COLORIMETRIC	U/L	115
MAGNESIUM	CALMAGITE	mEq/L	3.56
PHOSPHORUS	MOLYBDATE U. V	mg/dl	5.97
TOTAL PROTEIN	BIURET	g/dl	7.04
TRIGLYCERIDES	GPO/PAP	mg/dl	160
UREA	GLDH KINETIC	mg/dl	102
URIC ACID	URICASE /PAP	mg/dl	5.59

FOR IN VITRO DIAGNOSTIC USE ON

Manufactured By:
Coral Clinical Systems
A Division of Tulip Diagnostics (P) Ltd.
Building E, Plot No. M-46/47, Phase III B,
Verna Industrial Estate, Verna, Goa – 403 722, India.

REGD. OFFICE:
Gitanjali, Tulip Block,
Dr. Antonio Do Rego Bagh,
Alto Santacruz, Bambolim Complex P.O.,
Goa – 403 202, India.

ChemCAL™

Chemistry Calibrator

LOT 13925051A **IVD** **Mfd.** May. 2025 **Exp.** Mar. 2027

CORALYZER 100/ 200/ MINI

PRODUCT	METHOD	UNITS	CONCENTRATION
ALBUMIN	BROMOCRESOL GREEN	g/dl	4.76
ALKALINE PHOSPHATASE	IFCC (AMP)	U/L	413
	DGKC (DEA)	U/L	443
ALT/ GPT	MOD. IFCC	U/L	120
AST/ GOT	MOD. IFCC	U/L	121
AMYLASE	CNPG3, KINETIC	U/L	214
BILIRUBIN DIRECT	MOD. JENDRASSIK-GROF	mg/dl	1.93
	DCA		2.73
BILIRUBIN TOTAL	MOD. JENDRASSIK-GROF	mg/dl	6.79
	DCA		4.37
CALCIUM	OCPC	mg/dl	10.7
	AIII	mg/dl	10.9
CHLORIDE	THIOCYANATE	mmol/l	110
CHOLESTEROL	CHOD-PAP	mg/dl	218
HDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	65
LDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	85
CREATINE KINASE (CK)	MOD. IFCC	U/L	249
CREATININE	MOD. JAFFE'S KINETIC	mg/dl	3.38
	ENZYMATIC PAP	mg/dl	4.64
ɣ- GLUTAMYL TRANSFERASE (GGT)	CARBOXY SUBSTRATE	U/L	134
GLUCOSE	GOD/POD	mg/dl	192
IRON	FERROZINE	µg/dl	140
LACTATE DEHYDROGENASE	DGKC (P-L)	U/L	494
LIPASE	COLORIMETRIC	U/L	118
MAGNESIUM	CALMAGITE	mEq/L	3.36
PHOSPHORUS	MOLYBDATE U. V	mg/dl	5.57
TOTAL PROTEIN	BIURET	g/dl	7.84
TRIGLYCERIDES	GPO/PAP	mg/dl	173
UREA	GLDH KINETIC	mg/dl	119
URIC ACID	URICASE /PAP	mg/dl	5.84

FOR IN VITRO DIAGNOSTIC USE ONLY

Manufactured By:
Coral Clinical Systems
A Division of Tulip Diagnostics (P) Ltd.
Building E, Plot No. M-46/47, Phase III B,
Verna Industrial Estate, Verna, Goa – 403 722, India.

REGD. OFFICE:
Gitanjali, Tulip Block,
Dr. Antonio Do Rego Bagh,
Alto Santacruz, Bambolim Complex P.O.,
Goa – 403 202, India.