ChemCAL™

Chemistry Calibrator

LOT 13925122 **IVD** **Mfd.** DEC. 2025 **Exp.** JUL. 2027

BT 1500 / BT 4500

PRODUCT	METHOD	UNITS	CONCENTRATION
ALBUMIN	BROMOCRESOL GREEN	g/dl	4.93
ALKALINE PHOSPHATASE	IFCC (AMP)	U/L	364
ALT/ GPT	MOD. IFCC	U/L	95.0
AST/ GOT	MOD. IFCC	U/L	109
AMYLASE	CNPG3, KINETIC	U/L	199
BILIRUBIN DIRECT	MOD. JENDRASSIK - GROF	mg/dl	1.69
	DCA	mg/dl	2.00
BILIRUBIN TOTAL	MOD. JENDRASSIK - GROF	mg/dl	4.55
	DCA	mg/dl	4.20
CALCIUM	AIII	mg/dl	11.9
CHOLESTEROL	CHOD-PAP	mg/dl	212
HDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	54.0
LDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	141
CREATINE KINASE (CK)	MOD. IFCC	U/L	243
CK-MB	MOD. IFCC	U/L	98.0
CREATININE	MOD. JAFFE'S KINETIC	mg/dl	3.74
	ENZYMATIC PAP	mg/dl	4.59
γ- GLUTAMYL TRANSFERASE (GGT)	CARBOXY SUBSTRATE	U/L	130
GLUCOSE	GOD/POD	mg/dl	201
LACTATE DEHYDROGENASE	DGKC (P -L)	U/L	250
LIPASE	COLORIMETRIC	U/L	75.0
MAGNESIUM	CALMAGITE	mEq/L	3.00
PHOSPHORUS	MOLYBDATE U.V	mg/dl	5.69
TOTAL PROTEIN	BIURET	g/dl	7.99
TRIGLYCERIDES	GPO/PAP	mg/dl	181
UREA	GLDH KINETIC	mg/dl	96.0
URIC ACID	URICASE /PAP	mg/dl	6.33

FOR IN VITRO DIAGNOSTIC USE ONLY

Manufactured By: Coral Clinical Systems A Division of Tulip Diagnostics (P) Ltd. Building E, Plot No. M-46/47, Phase III B, Verna Industrial Estate, Verna, Goa – 403 722, India.

REGD. OFFICE:
Gitanjali, Tulip Block,
Dr. Antonio Do Rego Bagh,
Alto Santacruz, Bambolim Complex P.O.,
Goa – 403 202, India.

ChemCAL™

Chemistry Calibrator

LOT 13925122 **IVD** **Mfd.** DEC. 2025 **Exp.** JUL. 2027

CORALYZER SMART / JUNIOR

PRODUCT	METHOD	UNITS	CONCENTRATION
ALBUMIN	BROMOCRESOL GREEN	g/dl	4.80
ALKALINE PHOSPHATASE	IFCC (AMP)	U/L	357
ALT/ GPT	MOD. IFCC	U/L	99.0
AST/ GOT	MOD. IFCC	U/L	94.0
AMYLASE	CNPG3, KINETIC	U/L	197
BILIRUBIN DIRECT	DCA	mg/dl	2.60
BILIRUBIN TOTAL	DCA	mg/dl	4.50
CALCIUM	AIII	mg/dl	10.8
CHOLESTEROL	CHOD-PAP	mg/dl	205
HDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	63.0
LDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	121
CREATINE KINASE (CK)	MOD. IFCC	U/L	237
CK-MB	MOD. IFCC	U/L	91.0
CREATININE	MOD. JAFFE'S KINETIC	mg/dl	4.70
	ENZYMATIC PAP	mg/dl	4.40
γ- GLUTAMYL TRANSFERASE (GGT)	CARBOXY SUBSTRATE	U/L	122
GLUCOSE	GOD/POD	mg/dl	193
LACTATE DEHYDROGENASE	DGKC (P -L)	U/L	250
LIPASE	COLORIMETRIC	U/L	75.0
MAGNESIUM	CALMAGITE	mEq/L	3.20
PHOSPHORUS	MOLYBDATE U.V	mg/dl	5.60
TOTAL PROTEIN	BIURET	g/dl	7.40
TRIGLYCERIDES	GPO/PAP	mg/dl	165
UREA	GLDH KINETIC	mg/dl	105
URIC ACID	URICASE /PAP	mg/dl	6.00

FOR IN VITRO DIAGNOSTIC USE ONLY

Manufactured By: Coral Clinical Systems A Division of Tulip Diagnostics (P) Ltd. Building E, Plot No. M-46/47, Phase III B, Verna Industrial Estate, Verna, Goa – 403 722, India.

REGD. OFFICE: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa – 403 202, India.