Coral Clinical Systems
A Division of Tulip Diagnostics (P) Ltd.
Verna, Goa – 403 722, INDIA

Safety Data Sheet

1. INFORMATION ON THE SUBSTANCE / PREPARATION AND COMPANY

1.1 Product Name : **Acid Phosphatase Kit (Mod. King's method)**

Kit	Catalog Number	Components
Acid Phosphatase Kit (Mod. King's method)	1101010010	L1-Buffer Reagent L2-Substrate Reagent L3-Colour Reagent L4-Tartrate Reagent S-Phenol Standard

1.2 Intended Use : For the determination of Acid Phosphatase activity in serum.

1.3 Company : **Coral Clinical Systems**

(A Division of Tulip Diagnostics (Pvt) Ltd.)
Plot No. M-46 & M-47, Building E,
Phase III B, Verna Industrial Estate,
Verna-Goa, 403722
Tel :- +91-832-6680121
Fax :- +91-832-2887028
E-mail :- coral@tulipgroup.com

1.4 In emergencies Call your local emergency center.

2.Composition / Information on Hazardous Ingredients

Chemical Name	CAS #	% W/V	Exposure Limits in Air				
			ACGIH		OSHA		OTHER
			TLV	STEL	PEL	STEL	
Reagent 1: Buffer Reagent							
Sodium Hydroxide	1310-73-2	0.75 %	N / A	2 mg/m^{3}C	2 mg/m^3	N / A	NIOSH ILDH 10mg/m^3
Reagent 2: Substrate Reagent							
NA							
Reagent 3: Colour Reagent							
Sodium Hydroxide	1310-73-2	0.8 %	N / A	2 mg/m^{3}C	2 mg/m^3	N / A	NIOSH ILDH 10mg/m^3
Reagent 4: Tartrate Reagent							
NA							

Standard

NA							

3.Hazard Identification

According to 1999/45/EC, the preparation is classified as hazardous substance or mixture.
Skin Irritation,category 2; H315
Eye Irritation,Category 2;H319

Other Information
The reagents do not contain materials from animal origin

4.First Aid Measures

Inhalation:
If breathing becomes difficult, remove victim to fresh air. Seek medical attention immediately.

Ingestion:
Get medical attention if there has been ingestion of this product.

Skin Contact:
Avoid skin contact. If skin contact occurs, remove contaminated clothing and wash exposed skin with water. Get medical attention.

Eye Contact:
Immediately flush eye(s) with large volume of water for atleast 15 minutes, occasionally lifting the lower lids. Get medical attention immediately.

5.Fire Fighting Measures

Flash Point (Method used): NA **Flammable Limits – LEL:** NA **UEL:** NA

Extinguishing Media:
Use fire-extinguishing media appropriate for site conditions.

Special Fire Procedures:
Structural firefighting gear and self-contained breathing apparatus will provide adequate protection if this product is in a fire area.

Unusual Fire and Explosion Hazards:
Sodium azide can react with copper, lead, brass, or solder in plumbing to form explosive compounds of lead azide and copper azide. Sodium azide can react with acids to form explosive hydrogen azide.

6.Accidental Release Measures

Steps to be taken in case material is Released or Spilled:
Initiate cleaning up the spill with plenty of water and absorb material to contain / pick up the spilled solution. Place all contaminated disposable into a suitable container, seal, label and hold for disposal.

7.Handling and Storage

See package insert for storage information.

8.Exposure Controls and Personal Protection

Ventilation Data:
Supplemental ventilation such as general exhaust ventilation should be adequate under normal use of this product.

Respiratory Protection:
Respiratory protection is not required under normal use of this product. If respiratory protection is needed, follow OSHA respirator regulations (29CFR1910.134) and, if necessary, wear a NIOSH approved respirator. Select respirator based on its suitability to provide worker protection for given working conditions, level of airborne concentration, and presence of sufficient oxygen.

Protective Gloves:
Wear appropriate protective gloves to prevent skin contact. Replace torn or punctured gloves promptly.

Other Protective Equipment:
Wear appropriate eye protection to prevent eye contact. Wear appropriate body protection to prevent skin contact.

Other Engineering Controls:
Eye wash stations and deluge showers.

Work Practices:
Good laboratory technique should be used when handling this product. Observe appropriate chemical hygiene. Avoid contact with eyes or skin. Do not place in mouth.

Hygienic Practices:
Do not eat, drink, or smoke while working with product. Upon completion of work activities involving this product, wash hands thoroughly with soap and water.

9.Physical And Chemical Properties

L1 - Buffer Reagent

Relative Vapour density(air = 1)	:	NA	Evaporation rate(nBuAc = 1):		NA
Specific Gravity (water = 1)	:	NA	Freezing / Melting Point	:	NA
Solubility in Water	:	Soluble	Boiling Point	:	NA
Vapour Pressure, mm Hg @ 20^0C:		NA	pH	:	4.8±0.05

L2 - Substrate Reagent

Relative Vapour density(air = 1)	:	NA	Evaporation rate(nBuAc = 1):		NA
Specific Gravity (water = 1)	:	NA	Freezing / Melting Point :		NA
Solubility in Water	:	Soluble	Boiling Point	:	NA
Vapour Pressure, mm Hg @ 20^0C:		NA	pH	:	10.0±0.5

L3- Colour Reagent

Relative Vapour density(air = 1)	:	NA	Evaporation rate(nBuAc = 1):		NA
Specific Gravity (water = 1)	:	NA	Freezing / Melting Point :		NA
Solubility in Water	:	Soluble	Boiling Point	:	NA
Vapour Pressure, mm Hg @ 20^0C:		NA	pH	:	10.0±0.05

L4- Tartrate Reagent

Relative Vapour density(air = 1)	:	NA	Evaporation rate(nBuAc = 1):		NA
Specific Gravity (water = 1)	:	NA	Freezing / Melting Point :		NA
Solubility in Water	:	Soluble	Boiling Point	:	NA
Vapour Pressure, mm Hg @ 20^0C:		NA	pH	:	5.00±0.05

S- Phenol Standard

Relative Vapour density(air = 1)	:	NA	Evaporation rate(nBuAc = 1):		NA
Specific Gravity (water = 1)	:	NA	Freezing / Melting Point	:	NA
Solubility in Water	:	Soluble	Boiling Point	:	NA
Vapour Pressure, mm Hg @ 20^0C:		NA	pH	:	NA

Odour and Appearance Information
Reagent 1: Clear, colourless solution

Reagent 2: Clear to pale yellow solution

Reagent 3: Clear, yellow solution

Reagent 4: Clear, colourless solution

Reagent 5: Clear, colourless solution

10.Stability and Reactivity

Incompatibility (Materials to Avoid):
Strong bases, strong acids, and water reactive materials.

Hazardous Decomposition Products:
Thermal decomposition may produce carbon monoxide and carbon dioxide.

Will Hazardous Polymerization Occur?
Hazardous polymerization will not occur.

Conditions to Avoid / Polymerization: NA.

Is the Product Stable?
Yes, under normal handling and storage conditions.

Conditions to Avoid/stability
Stable Solution. Avoid acidification of solution, which may generate hydrogen cyanide gas.

11.Toxicological Information

Toxicity Data:
Sodium azide (undiluted): LD50 (rat and mouse, oral)=27mg/kg; LD50 (rabbit, skin)=20 mg/kg.

Reproductive effects:
NA.

Target organ Effects:
Eyes (redness), Skin (redness), central nervous systems (nausea/vomiting), cardiovascular systems (fall in blood pressure, change in heart rate), digestive (nausea/vomiting, diarrhea).

Carcinogenicity: No

CHEMICAL NAME	CAS #	% W/V	NTP Carcinogen		IARC	OSHA
			Known	Anticipated		
NA.						

12.Ecological Information

Environmental Fate / Stability: NA

Effect of Material on plants or animals: NA

Effect of Chemical on Aquatic Life: NA

13.Disposal Considerations

EPA Waste Number and Proper Waste Disposal Method:
Please consult local, state and federal regulations for additional guidance on disposal.

14.Transportation Information

Is this Material Hazardous? Not regulated under transportation regulations.

Proper Shipping Name : NA		
Hazard Class Number : NA	Packing Group: NA	UN Number: NA

15.Regulatory Information: NA.

16.Other Information

Hazard Statement and Precautionary statements referred elsewhere in the safety data sheet:
H315 - Causes skin irritation.
H319-Causes serious eye irritation.

The information provided on this MSDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text. It remains the user's own responsibility to make sure that the information is appropriate and complete for his specific use of this product. The user is also responsible for observing any laws and applicable guidelines.

Effective Date:02/09/2021 **Version No: 00**