Coral Clinical Systems
A Division of Tulip Diagnostics (P) Ltd.
Verna, Goa – 403 722, INDIA

Safety Data Sheet

1. INFORMATION ON THE SUBSTANCE / PREPARATION AND COMPANY

1.1 Product Name : **Calcium Kit (OCPC Method)**

Kit	Catalog No.	Components
Calcium Kit (OCPC Method)	1103010235/1103010275/1103012150/ 1103010010/1105050060/1105050240/ 1105050040/1106060240/1106060040/ 1116290250/1116290450/1123070060	L1-Buffer Reagent L2-Colour Reagent S-Calcium Standard

1.2 Intended Use : For the determination of Calcium in serum or plasma.

1.3 Company : **Coral Clinical Systems**
(A Division of Tulip Diagnostics (Pvt) Ltd.)
Plot No. M-46 & M-47, Building E,
Phase III B, Verna Industrial Estate,
Verna-Goa, 403722
Tel :- +91-832-6680121
Fax :- +91-832-2887028
E-mail :- coral@tulipgroup.com

1.4 In emergencies :Call your local emergency center.

2.Composition / Information on Hazardous Ingredients

Chemical Name	CAS #	% W/V	Exposure Limits in Air				
			ACGIH		OSHA		OTHER
			TLV	STEL	PEL	STEL	

Reagent 1: **No Hazardous components**

Reagent 2: **No Hazardous components**

Reagent 3: **No Hazardous components**

3.Hazard Identification

According to 1999/45/EC, the preparation is classified as non hazardous substance or mixture.

Other Information
The reagents do not contain materials from animal origin

4.First Aid Measures

Inhalation:
If breathing becomes difficult, remove victim to fresh air. Seek medical attention immediately.

Ingestion:
Get medical attention if there has been ingestion of this product.

Skin Contact:
Avoid skin contact. If skin contact occurs, remove contaminated clothing and wash exposed skin with water for atleast 15 minutes. Get medical attention immediately.

Eye Contact:
Immediately flush eye(s) with large volume of water for atleast 15 minutes, occasionally lifting the lower lids. Get medical attention immediately.

5.Fire Fighting Measures

Flash Point (Method used): NA **Flammable Limits – LEL:** NA **UEL:** NA

Extinguishing Media:
Use fire extinguishing media appropriate for site conditions.

Special Fire Procedures:
Structural firefighting gear and self-contained breathing apparatus will provide adequate protection if this product is in a fire area.

Unusual Fire and Explosion Hazards:
Sodium azide can react with copper, lead, brass, or solder in plumbing to form explosive compounds of lead azide and copper azide. Sodium azide can react with acids to form explosive hydrogen azide.

6.Accidental Release Measures
Steps to be taken in case material is Released or Spilled:
PPE should be level D: lab gloves, chemical resistant apron, boots and splash goggles. Use an absorbent material to contain / pick up the spilled solution. Place all spill residues into a suitable container, seal, label and hold for disposal.

7.Handling and Storage
Refer to packet insert for additional information on handling and storage procedures.

8.Exposure Controls and Personal Protection

Ventilation Data:
A system of local and / or general exhaust is recommended to keep employee exposures as low as possible. Local exhaust ventilation is generally preferred because it can control the emissions of the contaminant at its source.

Respiratory Protection:
Respiratory protection is not required under normal use of this product. If respiratory protection is needed, follow OSHA respirator regulations (29CFR1910.134) and, if necessary, wear a NIOSH approved respirator. Select respirator based on its suitability to provide worker protection for given working conditions, level of airborne concentration, and presence of sufficient oxygen.

Protective Gloves:
Wear appropriate protective gloves to prevent skin contact. Replace torn or punctured gloves promptly.

Other Protective Equipment:
Wear appropriate eye protection to prevent eye contact. Wear appropriate body protection to prevent skin contact.

Other Engineering Controls:
Eye wash stations and deluge showers.

Work Practices:
Good laboratory technique should be used when handling this product. Observe appropriate chemical hygiene. Avoid contact with eyes or skin. Do not place in mouth.

Hygienic Practices:
Do not eat, drink, or smoke while working with product. Upon completion of work activities involving this product, wash hands thoroughly with soap and water.

9.Physical And Chemical Properties

L1- Buffer Reagent

Relative Vapour density(air = 1)	: NA	Evaporation rate(nBuAc = 1):	NA
Specific Gravity (water = 1)	: NA	Freezing / Melting Point	: NA
Solubility in Water	: Soluble	Boiling Point	: NA
Vapour Pressure, mm Hg @ 200C:	NA	pH	: NA

L2 - Colour Reagent

Relative Vapour density(air = 1)	: NA	Evaporation rate(nBuAc = 1):	NA
Specific Gravity (water = 1)	: NA	Freezing / Melting Point	: NA
Solubility in Water	: Soluble	Boiling Point	: NA
Vapour Pressure, mm Hg @ 200C	: NA	pH	: NA

S - Calcium Standard

Relative Vapour density(air = 1)	: NA	Evaporation rate(nBuAc = 1) :	NA
Specific Gravity (water = 1)	:NA	Freezing / Melting Point	: NA
Solubility in Water	:Soluble	Boiling Point	: NA
Vapour Pressure, mm Hg @ 200C	:NA	pH	: NA

Odour and Appearance Information
Reagent 1: Clear, colorless solution

Reagent 2: Clear, yellow solution

Reagent 3: Clear, colorless solution

10. Stability and Reactivity

Incompatibility (Materials to Avoid):
Strong bases, strong acids, and water reactive materials.

Hazardous Decomposition Products:
Thermal decomposition may produce carbon monoxide and carbon dioxide.

Will Hazardous Polymerization Occur?
Hazardous polymerization will not occur.

Conditions to Avoid / Polymerization: NA

Is the Product Stable?
Yes, under normal handling and storage conditions.

Conditions to Avoid/stability
Stable Solution. Avoid acidification of solution, which may generate hydrogen cyanide gas.

11.Toxicological Information

Toxicity Data:
Sodium Azide(Undiluted): LD50(rat and mouse, oral)=27mg/kg; LD50(rabbit, skin)=20mg/kg.
Hydrochloric acid: eye irritancy (rabbit)-1% solution/20sec., comeal scarring; Eye irritancy (rabbit)= 5mg/30sec, mild; eye irritancy (rabbit)=100mg/rinse, mild; skin irritancy(rabbit)=0.5ml/17% solution/4hours, corrosive burns.

Reproductive effects:
NA.

Target organ Effects:
Sodium Azide: Eyes, Skin, central nervous systems, cardiovascular systems. **Hydrochloric Acid**: Eyes, Skin and Respiratory system. Sodium Hydroxide: Eyes, Skin and Respiratory system.

Carcinogenicity: No

CHEMICAL NAME	CAS #	% W/V	NTP Carcinogen		IARC	OSHA
			Known	Anticipated		
NA.						

12.Ecological Information

Environmental Fate / Stability :NA

Effect of Material on plants or animals :NA

Effect of Chemical on Aquatic Life :NA

13..Disposal Considerations

EPA Waste Number and Proper Waste Disposal Method:
Please consult local, state and federal regulations for additional guidance on disposal.

14.Transportation Information

Is this Material Hazardous? Not regulated under transportation regulations.

Proper Shipping Name : NA		
Hazard Class Number : NA	Packing Group: NA	UN Number: NA

15.Regulatory Information NA.

16.Other Information

This product is designed for use by professionals.

The information provided on this MSDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text. It remains the user's own responsibility to make sure that the information is appropriate and complete for his specific use of this product. The user is also responsible for observing any laws and applicable guidelines.

Effective Date:06/01/2023 **Version No: 00**