ChemCAL™

Chemistry Calibrator

LOT 13926072 **IVD** **Mfd.** JUL. 2026 **Exp.** JAN. 2028

BT 1500 / BT 4500

PRODUCT	METHOD	UNITS	CONCENTRATION
ALBUMIN	BROMOCRESOL GREEN	g/dl	3.9
ALKALINE PHOSPHATASE	IFCC (AMP)	U/L	244
ALT/ GPT	MOD. IFCC	U/L	96
AST/ GOT	MOD. IFCC	U/L	106
AMYLASE	CNPG3, KINETIC	U/L	207
BILIRUBIN DIRECT	MOD. JENDRASSIK GROF	mg/dl	2.6
	DCA	mg/dl	2.4
BILIRUBIN TOTAL	MOD. JENDRASSIK GROF	mg/dl	4.7
	DCA	mg/dl	4.5
CALCIUM	AIII	mg/dl	10.5
CHOLESTEROL	CHOD-PAP	mg/dl	173
CHLORIDE	THIOCYNATE	Mmol/l	104
HDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	58
LDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	68
CREATINE KINASE (CK)	MOD. IFCC	U/L	348
CK-MB	MOD. IFCC	U/L	371
CREATININE	MOD. JAFFE'S KINETIC	mg/dl	3.5
	ENZYMATIC PAP	mg/dl	3.9
ɣ- GLUTAMYL TRANSFERASE (GGT)	CARBOXY SUBSTRATE	U/L	82
GLUCOSE	GOD/POD	mg/dl	204
LACTATE DEHYDROGENASE	DGKC (P -L)	U/L	548
LIPASE	COLORIMETRIC	U/L	101
MAGNESIUM	CALMAGITE	mEq/L	2.7
PHOSPHORUS	MOLYBDATE U.V	mg/dl	5.7
TOTAL PROTEIN	BIURET	g/dl	5.5
TRIGLYCERIDES	GPO/PAP	mg/dl	151
UREA	GLDH KINETIC	mg/dl	102
URIC ACID	URICASE /PAP	mg/dl	5.5
IRON	FERROZINE	µg/dl	84

FOR IN VITRO DIAGNOSTIC USE ONLY

Manufactured By: Coral Clinical Systems A Division of Tulip Diagnostics (P) Ltd. Building E, Plot No. M-46/47, Phase III B, Verna Industrial Estate, Verna, Goa – 403 722, India.

REGD. OFFICE:
Gitanjali, Tulip Block**,**
Dr. Antonio Do Rego Bagh,
Alto Santacruz, Bambolim Complex P.O.,
Goa – 403 202, India

ChemCAL™

Chemistry Calibrator

LOT 13926072 **IVD** **Mfd.** JUL. 2026 **Exp.** JAN. 2028

CORALYZER SMART / JUNIOR

PRODUCT	METHOD	UNITS	CONCENTRATION
ALBUMIN	BROMOCRESOL GREEN	g/dl	3.8
ALKALINE PHOSPHATASE	IFCC (AMP)	U/L	246
ALT/ GPT	MOD. IFCC	U/L	102
AST/ GOT	MOD. IFCC	U/L	103
AMYLASE	CNPG3, KINETIC	U/L	200
BILIRUBIN DIRECT	DCA	mg/dl	2.7
BILIRUBIN TOTAL	DCA	mg/dl	4.9
CALCIUM	AIII	mg/dl	10.1
CHOLESTEROL	CHOD-PAP	mg/dl	167
CHLORIDE	THIOCYNATE	Mmol/l	84
HDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	55
LDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	95
CREATINE KINASE (CK)	MOD. IFCC	U/L	385
CK-MB	MOD. IFCC	U/L	304
CREATININE	MOD. JAFFE'S KINETIC	mg/dl	3.3
	ENZYMATIC PAP	mg/dl	3.5
ɣ- GLUTAMYL TRANSFERASE (GGT)	CARBOXY SUBSTRATE	U/L	87
GLUCOSE	GOD/POD	mg/dl	196
LACTATE DEHYDROGENASE	DGKC (P-L)	U/L	556
LIPASE	COLORIMETRIC	U/L	103
MAGNESIUM	CALMAGITE	mEq/L	2.4
PHOSPHORUS	MOLYBDATE U.V	mg/dl	5.1
TOTAL PROTEIN	BIURET	g/dl	5.1
TRIGLYCERIDES	GPO/PAP	mg/dl	148
UREA	GLDH KINETIC	mg/dl	99
URIC ACID	URICASE /PAP	mg/dl	5.1
IRON	FERROZINE	µg/dl	54

FOR IN VITRO DIAGNOSTICS USE ONLY

Manufactured By: Coral Clinical Systems A Division of Tulip Diagnostics (P) Ltd. Building E, Plot No. M-46/47, Phase III B, Verna Industrial Estate, Verna, Goa – 403 722, India.

REGD. OFFICE:
Gitanjali, Tulip Block,
Dr. Antonio Do Rego Bagh,
Alto Santacruz, Bambolim Complex P.O.,
Goa – 403 202, India