ChemCAL™

Chemistry Calibrator

LOT 13926093 **IVD** **Mfd.** SEP 2026 **Exp.** FEB.2028

BT 1500 / BT 4500

PRODUCT	METHOD	UNITS	CONCENTRATION
ALBUMIN	BROMOCRESOL GREEN	g/dl	4.0
ALKALINE PHOSPHATASE	IFCC (AMP)	U/L	254
ALT/ GPT	MOD. IFCC	U/L	105
AST/ GOT	MOD. IFCC	U/L	95
AMYLASE	CNPG3, KINETIC	U/L	176
BILIRUBIN DIRECT	MOD. JENDRASSIK GROF	mg/dl	2.0
	DCA	mg/dl	1.8
BILIRUBIN TOTAL	MOD. JENDRASSIK GROF	mg/dl	4.7
	DCA	mg/dl	4.8
CALCIUM	AIII	mg/dl	11
CHOLESTEROL	CHOD-PAP	mg/dl	158
CHLORIDE	THIOCYNATE	Mmol/l	82
HDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	54
LDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	-
CREATINE KINASE (CK)	MOD. IFCC	U/L	344
CK-MB	MOD. IFCC	U/L	-
CREATININE	MOD. JAFFE'S KINETIC	mg/dl	3.7
	ENZYMATIC PAP	mg/dl	3.8
ɣ- GLUTAMYL TRANSFERASE (GGT)	CARBOXY SUBSTRATE	U/L	90
GLUCOSE	GOD/POD	mg/dl	198
LACTATE DEHYDROGENASE	DGKC (P -L)	U/L	428
LIPASE	COLORIMETRIC	U/L	108
MAGNESIUM	CALMAGITE	mEq/L	2.4
PHOSPHORUS	MOLYBDATE U.V	mg/dl	5.1
TOTAL PROTEIN	BIURET	g/dl	5.3
TRIGLYCERIDES	GPO/PAP	mg/dl	145
UREA	GLDH KINETIC	mg/dl	101
URIC ACID	URICASE /PAP	mg/dl	5.2
IRON	FERROZINE	µg/dl	196

FOR IN VITRO DIAGNOSTIC USE ONLY

Manufactured By: Coral Clinical Systems A Division of Tulip Diagnostics (P) Ltd. Building E, Plot No. M-46/47, Phase III B, Verna Industrial Estate, Verna, Goa – 403 722, India.

REGD. OFFICE:
Gitanjali, Tulip Block,
Dr. Antonio Do Rego Bagh,
Alto Santacruz, Bambolim Complex P.O.,
Goa – 403 202, India

ChemCAL™

Chemistry Calibrator

LOT 13926093 **IVD** **Mfd.** SEP. 2026 **Exp.** FEB. 2028

CORALYZER SMART / JUNIOR

PRODUCT	METHOD	UNITS	CONCENTRATION
ALBUMIN	BROMOCRESOL GREEN	g/dl	3.9
ALKALINE PHOSPHATASE	IFCC (AMP)	U/L	253
ALT/ GPT	MOD. IFCC	U/L	107
AST/ GOT	MOD. IFCC	U/L	96
AMYLASE	CNPG3, KINETIC	U/L	177
BILIRUBIN DIRECT	DCA	mg/dl	2.5
BILIRUBIN TOTAL	DCA	mg/dl	4.7
CALCIUM	AIII	mg/dl	11
CHOLESTEROL	CHOD-PAP	mg/dl	158
CHLORIDE	THIOCYNATE	Mmol/l	85
HDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	53
LDL-CHOLESTEROL	DIRECT ENZYMATIC	mg/dl	-
CREATINE KINASE (CK)	MOD. IFCC	U/L	342
CK-MB	MOD. IFCC	U/L	-
CREATININE	MOD. JAFFE'S KINETIC	mg/dl	3.5
	ENZYMATIC PAP	mg/dl	3.8
ɣ- GLUTAMYL TRANSFERASE (GGT)	CARBOXY SUBSTRATE	U/L	87
GLUCOSE	GOD/POD	mg/dl	197
LACTATE DEHYDROGENASE	DGKC (P-L)	U/L	420
LIPASE	COLORIMETRIC	U/L	109
MAGNESIUM	CALMAGITE	mEq/L	2.2
PHOSPHORUS	MOLYBDATE U.V	mg/dl	5.1
TOTAL PROTEIN	BIURET	g/dl	5.2
TRIGLYCERIDES	GPO/PAP	mg/dl	140
UREA	GLDH KINETIC	mg/dl	102
URIC ACID	URICASE /PAP	mg/dl	5.2
IRON	FERROZINE	µg/dl	196

FOR IN VITRO DIAGNOSTICS USE ONLY

Manufactured By: Coral Clinical Systems A Division of Tulip Diagnostics (P) Ltd. Building E, Plot No. M-46/47, Phase III B, Verna Industrial Estate, Verna, Goa – 403 722, India.

REGD. OFFICE:
Gitanjali, Tulip Block,
Dr. Antonio Do Rego Bagh,
Alto Santacruz, Bambolim Complex P.O.,
Goa – 403 202, India